USER MANUAL

Placing the Essence Air™ mattress overlay in use

- Remove the Essence Air™ out of the packaging. Do not use knives, scissors or other sharp objects to open the package.
- Remove all bed linen.
- Place the Essence Air™ on the base mattress at the foot of the bed.
- Open both valves and unroll the Essence Air™ until it is flat on the base mattress.

To inflate the Essence Air™

Turn both valves at the foot end of the Essence Air™ in anti-clockwise direction. The Essence Air™ will self-inflate in approximately 3-5 minutes. When the Essence Air™ has inflated, close both valves firmly to achieve a full seal.

Note

Full self-inflation on first use requires minimum 10 minutes with valves fully opened.

Deflating after use

- Open both valves completely
- Start rolling from the non-valve end.
- Roll by hand as tightly as possible in the direction of the valve end of the Essence Air™.
- You will hear the air exhausting through the valves.
- When the Essence Air™ is fully rolled and compressed, close the valves securely.
- The Essence Air™ is ready for transport.

Cleaning and care

- Hygienic wipes can be used or the Essence Air™ can be wiped down with water and any mild detergent. Please ensure the correct dilution is used.
- Never use phenol based cleaning agents or organic solvents to clean the Essence Air™.
- Do not tumble dry.
- Do not dry clean.
- Do not machine wash.

Caution

- Never apply direct flame/heat to the Essence Air™.
- Do not touch the Essence Air™ with any sharp object.

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ALL RIGHTS RESERVED